

Document Number:	TL19	Version Number:	3.0
Written by:	WADA Science	Approved by:	WADA Executive Committee
Reviewed by:	WADA Laboratory Expert Group		
Date:	21 December 2020	Effective Date:	1 January 2021

IN SITU FORMATION OF PREDNISONE AND PREDNISOLONE

1.0 Introduction

WADA wishes to draw the attention of the Laboratories to the possible detection of **Prednisone** and/or **Prednisolone** in urine *Samples* resulting from *in situ* microbial transformation of **Cortisone** and **Cortisol**, respectively.

The endogenous glucocorticoids cortisol (also known as hydrocortisone) and cortisone have a 3-oxo-4-ene structure (Figure 1). This structural element is a possible microbial Δ^1 -steroid-dehydrogenase (Δ^1 -SDH) target for 1,2-dehydrogenation, which may transform cortisol and cortisone into prednisolone and prednisone, respectively [1].




Figure 1. The microbial conversion of cortisol to prednisolone by Δ^1 -SDH enzymatic activity.

2.0 Analysis and Reporting Requirements

Before reporting a result as an *Adverse Analytical Finding (AAF)* for prednisone and/or prednisolone, Laboratories shall evaluate whether the finding is the result of the *in situ* microbial transformation of cortisone and cortisol, respectively.

1. When prednisone and/or prednisolone are present at levels higher than (>) the *Minimum Reporting Limit (MRL)* of 30 ng/mL, perform a CP using an extraction method [e.g., Solid Phase Extraction (SPE)] prior to the enzymatic hydrolysis step in order to avoid inducing the *in situ* formation of prednisone and prednisolone by the enzymatic activity of microbes already present in the *Sample*.

[Comment: However, if the in situ formation of prednisone and prednisolone have already occurred prior to the enzymatic hydrolysis, SPE will have no impact.]

2. Evaluate the overall pattern of the *Metabolites* present in the *Sample*: the presence of aglycones only (*i.e.*, absence of the glucuronide *Metabolites*) and the absence of the 20 β -hydroxy-*Metabolite* can be used as an indication of microbial activity.
3. Perform GC/C/IRMS [2] analysis (depending on Laboratory's analytical capacity, which may require the subcontracting of the analysis) when prednisone and/or prednisolone are present at levels (after correction for SG, if needed) estimated between (\geq) 30 ng/mL and (\leq) 60 ng/mL.

WADA Technical Letter – TL19

Document Number:	TL19	Version Number:	3.0
Written by:	WADA Science	Approved by:	WADA Executive Committee
Reviewed by:	WADA <u>Laboratory Expert Group</u>		
Date:	21 December 2020	Effective Date:	1 January 2021

4. Report the result as a Negative Finding if:

- Prednisone and/or prednisolone are detected at an estimated concentration (SG-adjusted, if needed), which is lower than (<) the *MRL* of 30 ng/mL; or
- The estimated concentration of prednisone and/or prednisolone (SG-adjusted, if needed) is between (\geq) the *MRL* of 30 ng/mL and (\leq) 60 ng/mL, but:
 - The overall pattern of relevant substances and *Metabolites* present in the *Sample* indicates microbial activity (e.g., absence of 20 β -hydroxy- and glucuronide-*Metabolites* and/or indication of a Δ^1 -SDH activity); and/or
 - The GC/C/IRMS analysis demonstrates an endogenous origin of the substance(s) ^[2].

5. Report the result as an *Atypical Finding (ATF)* if:

- The estimated concentration of prednisone and/or prednisolone (SG-adjusted, if needed) is between (\geq) 30 ng/mL and (\leq) 60 ng/mL, and the GC/C/IRMS analysis is inconclusive or cannot be performed ^[2];
- The estimated concentration (SG-adjusted, if needed) of prednisone and/or prednisolone in the *Sample* is higher than (>) 60 ng/mL, but the *Sample* shows signs of extensive degradation;
- If an *ATF* is reported, the Laboratory shall include a comment in the *ADAMS* Test Report recommending the Testing Authority to conduct at least one (1) follow-up no-notice test on the *Athlete* within a reasonable time frame (e.g. within 2 weeks).

6. Report the result as an *Adverse Analytical Finding (AAF)* if:

- The confirmed estimated concentration (SG-adjusted, if needed) of prednisone and/or prednisolone in the *Sample* is higher than (>) 60 ng/mL, and there are no signs of extensive degradation; or
- The estimated concentration of prednisone and/or prednisolone (SG-adjusted, if needed) is between (\geq) 30 ng/mL and (\leq) 60 ng/mL, and the GC/C/IRMS analysis demonstrates an exogenous origin of the substance(s) ^[2];
- It is recommended that the Laboratory seeks a second opinion, in writing, from another Laboratory before reporting the *AAF*. The second opinion shall be recorded in the Laboratory Documentation Package.

3.0 References

[1] Bredehöft M., *et. al.* Investigations of the microbial transformation of cortisol to prednisolone in urine samples. *J Steroid Biochem Mol Biol.* **129**(1-2): 54-60, 2012.

[2] WADA *Technical Document* TD IRMS: Detection of Synthetic Forms of *Prohibited Substances* by GC/C/IRMS.

[Current versions of WADA *Technical Documents* may be found at <https://www.wada-ama.org/en/what-we-do/science-medical/laboratories>]