

## TD2019GH

### Summary of Major Modifications

The Technical Document on Human Growth Hormone (hGH) Isoform Differential Immunoassays for *Doping Control* Analyses, TD2019GH, has been modified to account for new analytical evidence showing that BD Vacutainer® SST™ tubes, US ref 367986 are Fit-for-Purpose<sup>1</sup>. Therefore, these tubes can be used for blood *Sample* collection for hGH analysis.

The TD2019GH replaces the former TD2015GH.

---

<sup>1</sup> Geoffrey D. Miller, H.D. Cox, V. Nair and D. Eichner. Anti-doping analytes in serum: a comparison of SST and SST-II Advance blood collection tubes. *Drug Test Anal* 2019 Mar 22 doi: 10.1002/dta.2592. [Epub ahead of print].